

Product datasheet for **RC207174A1V**

Human ULBP1 (NM_025218) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ULBP1 (NM_025218) AAV Particle
Tag:	Myc-DDK
Symbol:	ULBP1
Synonyms:	N2DL-1; NKG2DL1; RAET1I
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207174 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGCGGCCGCCAGCCCCGCTTCCTTCTGTGCCTCCCGCTTCTGCACCTGCTGTCTGGCTGGTCCC
GGGCAGGATGGGTCGACACACTGTCTTTGCTATGACTTCATCATCACTCCTAAGTCCAGACCTGAACC
ACAGTGGTGTGAAGTTCAAGGCCTGGTGGATGAAAGGCCTTTTCTTCACTATGACTGTGTTAACCACAAG
GCCAAAGCCTTTGCTTCTCTGGGAAGAAAGTCAATGTCACAAAACCTGGGAAGAACAACCTGAAACAC
TAAGAGACGTGGTGGATTTCTTAAAGGGCACTGCTTGACATTCAAGTGGAGAATTTAATACCCATTGA
GCCCCCTACCCTGCAGGCCAGGATGTCTTGAGCATGAAGCCCATGGACACGGCAGAGGATCTTGGCAG
TTCCTCTTCAATGGACAGAAGTTCCTCCTCTTTGACTCAAACAACAGAAAGTGGACAGCACTTCATCCTG
GAGCCAAGAAGATGACAGAGAAGTGGGAGAAGAACAGGGATGTGACCATGTTCTTCCAGAAGATTTCACT
GGGGGATTGTAAGATGTGGCTTGAAGAATTTTGTGACTGGGAACAAATGCTGGATCCAACAAAACCA
CCCTCTCTGGCCCCAGGCACAACCAACCAAGGCCATGGCCACCACCCTCAGTCCCTGGAGCCTTCTCA
TCATCTTCTCTGCTTCATTCTAGCTGGCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207174 protein sequence Red=Cloning site Green=Tags(s) MAAAASPALLCLPLLHLLSGWSRAGWVDTHCLCYDFIITPKSRPEPQWCEVQGLVDERPFLHYDCVNHK AKAFASLGKKVNVTKTWEEQTETLRDVDFLKGQLLDIQVENLPIEPLTLQARMSCEHAHGHGRGSWQ FLFNGQKFLLFDSNNRKWTALHPGAKKMTKEWKNRDVTMFFQKISLGDCKMWLEEFMYWEQMLDPTKP PSLAPGTTQPKAMATTLSPWSLLIIFLCFILAGR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_025218
ORF Size:	732 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_025218.2</u>
RefSeq Size:	3160 bp
RefSeq ORF:	735 bp
Locus ID:	80329
UniProt ID:	<u>Q9BZM6</u>
Cytogenetics:	6q25.1
MW:	28 kDa