

## Product datasheet for **RC207188L1V**

### DZIP1 (NM\_198968) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DZIP1 (NM_198968) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DZIP1
Synonyms:	DZIP; DZIPT1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_198968
ORF Size:	2601 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207188).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198968.2</a>
RefSeq Size:	7498 bp
RefSeq ORF:	2604 bp


[View online »](#)

Locus ID:	22873
UniProt ID:	<a href="#">Q86YF9</a>
Cytogenetics:	13q32.1
MW:	98.7 kDa
Gene Summary:	May participate in spermatogenesis via its interaction with DAZ (PubMed:15081113). Has a role in primary cilium formation (PubMed:19852954).[UniProtKB/Swiss-Prot Function]