

Product datasheet for **RC207189L4V**

RHOBTB3 (NM_014899) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RHOBTB3 (NM_014899) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RHOBTB3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014899
ORF Size:	1833 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207189).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014899.3 , NP_055714.2
RefSeq Size:	5566 bp
RefSeq ORF:	1836 bp
Locus ID:	22836
UniProt ID:	O94955
Cytogenetics:	5q15
Domains:	BTB
MW:	69.4 kDa


[View online »](#)

Gene Summary:

RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004]