

## Product datasheet for **RC207201L4V**

### NMUR1 (NM\_006056) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NMUR1 (NM_006056) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NMUR1
Synonyms:	(FM-3); FM-3; FM3; GPC-R; GPR66; NMU1R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006056
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207201).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006056.2</a>
RefSeq Size:	3274 bp
RefSeq ORF:	1281 bp
Locus ID:	10316
UniProt ID:	<a href="#">Q9HB89</a>
Cytogenetics:	2q37.1
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)

**MW:** 47.3 kDa

**Gene Summary:** Receptor for the neuromedin-U and neuromedin-S neuropeptides.[UniProtKB/Swiss-Prot Function]