

Product datasheet for **RC207217L2V**

PKC beta 1 (PRKCB) (NM_002738) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PKC beta 1 (PRKCB) (NM_002738) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PKC beta 1 |
| Synonyms: | PKC-beta; PKCB; PKCbeta; PKCI(2); PRKCB1; PRKCB2 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_002738 |
| ORF Size: | 2019 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC207217). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_002738.5 |
| RefSeq Size: | 8014 bp |
| RefSeq ORF: | 2022 bp |
| Locus ID: | 5579 |
| UniProt ID: | P05771 |
| Cytogenetics: | 16p12.2-p12.1 |
| Domains: | C2, pkinase, S_TK_X, TyrKc, DAG_PE-bind, S_TKc |
| Protein Families: | Druggable Genome, Protein Kinase |


[View online »](#)

| | |
|--------------------------|--|
| Protein Pathways: | B cell receptor signaling pathway, Calcium signaling pathway, Chemokine signaling pathway, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Gap junction, Glioma, GnRH signaling pathway, Leukocyte transendothelial migration, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Natural killer cell mediated cytotoxicity, Non-small cell lung cancer, Pathways in cancer, Phosphatidylinositol signaling system, Tight junction, Vascular smooth muscle contraction, VEGF signaling pathway, Vibrio cholerae infection, Wnt signaling pathway |
| MW: | 77 kDa |
| Gene Summary: | Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role in cells. The protein encoded by this gene is one of the PKC family members. This protein kinase has been reported to be involved in many different cellular functions, such as B cell activation, apoptosis induction, endothelial cell proliferation, and intestinal sugar absorption. Studies in mice also suggest that this kinase may also regulate neuronal functions and correlate fear-induced conflict behavior after stress. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008] |