

Product datasheet for **RC207234**

C19orf18 (NM_152474) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C19orf18 (NM_152474) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C19orf18
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC207234 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAAGGTTTCAGAGTGGTTTCCTCATTTTGTTTTTGTTTTAATGGAATGCCAACTTCATTTATGCT
TGCCGTATGCAGATGGACTCCATCCCACTGAAACATAACAGGCTTACCAGGCAGTAAACGGTCACAACC
ACCCAGAAACATCACAAAGAGCCCAAAGTGTTCTTTCTATAAAACCCAGTTGCCTGGGATTCAAGGGCT
GCCTCGAGATCCACGGCTGCATCCCCTACGAACCCCATGAAATTCCTGAGGAATAAAGCAATAATTCGGC
ATAGACCTGCTCTTGTTAAAGTAATTTAATTTTCGAGCGTAGCCCTCAGCATTGCCCTGATGTGGGAT
GGCAATCTCTATATGATATATCGACTGGCACAGGCTGAGGAAAGACAACAGCTCGAGTCATTTATAAG
AACCTCAGGATACCGTTATTAGGAGATGAAGAAGAGGGCTCAGAGGACGAGGGTGAGTCCACGCACCTAC
TTCCAGAGAACGAAAATGAGCTGAAAAAGTTTCATCCACTCAGTTATTATATCAAAAAGAAGCAAAAATAT
TAAGAAGAACTGAAGGAAGAGCAAACTCAGTAACAGAAAACAAAAGAATGCGTACATAATGGA
AAAATGGAAGACTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207234 protein sequence
Red=Cloning site Green=Tags(s)

MDKVQSGFLILFLMECQLHLCLPYADGLHPTGNITGLPGSKRSQPPRNITKEPKVFFHKTLPLGIQGA
ASRSTAASPTNPMKFLRNKAIIRHRPALVKVILISSVAFSIALICGMAISYMIYRLAQAEERQQLLESLYK
NLRIPLLGDEEEGSEDEGESTHLLPENENELEKFIHSVVIISKRSKNIKKKLKEEQNSVTENKTKNASHNG
KMEDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Chromatograms: https://cdn.origene.com/chromatograms/mk6341_a09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152474

ORF Size: 645 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152474.5](#)

RefSeq Size: 928 bp

RefSeq ORF: 648 bp

Locus ID: 147685


UniProt ID: [Q8NEA5](#)

Cytogenetics: 19q13.43


Protein Families: Transmembrane

MW: 24.2 kDa

Product images:



Circular map for RC207234



Western blot validation of overexpression lysate (Cat# [LY407528]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC207234 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).