

Product datasheet for **RC207235A1V**

Human ZNF785 (NM_152458) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZNF785 (NM_152458) AAV Particle
Tag:	Myc-DDK
Symbol:	ZNF785
Synonyms:	ZNF688
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207235 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGCCGCCCTGGCGCCGCGCCCGGCCACGTACCGGGGAGGCGGGCCCCGGAGGACGAGGGAAA
GCAGGCCGGCGCCGTGAGCTTCGCGGACGTGGCCGTGTAATTCTCTCCCGAGGAGTGGGAATGCCTGCG
GCCAGCGCAGAGGGCCCTGTACCGGGACGTGATGCGGGAGACCTTCGGCCACCTGGGCGCGCTGGGATTT
TCAGTTCCCAAACCCGCTTTTATCTCGTGGGTGGAAGGAGAAGTGGAGGCGTGGAGCCCGAGGCCCAGG
ATCCCGACGGTGAGAGCTCTGCAGCTTTCAGCAGGGGCCAAGGACAGGAAGCAGGATCCAGGGATGGGAA
TGAGGAGAAGGAAAGGCTGAAGAAGTGCCAAAACAAAAGAGGTGGCGCATGAAGTGGCTGTCAAGGAG
TGGTGGCCCGAGCGTCGCCTGCCAGAGTTCTGCAACCCTAGGCAGAGCCCATGAATCCCTGGCTCAAGG
ACACTCTGACCCGAAGACTGCCCCACTTTGCCAGACTGTGGCCGCACTTCAGCTACCCCTCCCTCCT
GGCCAGCCACCAGCGGTCCACTCCGGGAGCGGCCCTTCTCCTGCGGCCAGTGTGAGCGCGTTTCTCC
CAGCGCAGGTACCTGCTCCAGCATCAGTTTATCCACACCGGCGAGAAGCCCTACCCCTGCCCGACTGCG
GGCGCCGCTTCCGCCAGAGGGTTCCCTGGCTATCCACAGGCGGGCTCACACGGGGAGAAGCCTTACGC
GTGCTCAGACTGCAAGAGTCGTTCACTTACCCTACCTGCTGCGCATCCACAGCGCAAGCACACGGGC
GAGAAGCCCTACAGTGCCCCGATTGCAGCCTCCGTTTCGCTACACCTCCCTGCTGGCCATCCACAGGC
GCATACACACCGGCGAGAAGCCCTACCCCTGTCTGACTGCGGCGCGCGTTACCTATTCTTCCCTCCT
CCTCAGTCACCGGCGCATTCACTCCGACAGCCGCCCTTCCCCTGCGTGGAGTGTGGAAAGGCTTCAAG
CGCAAGACCGCCTGGAAGCCCATCGGTGGATCCACCGCTCCTGCAGCGAGAGGCGCGCTGGCAGCAGG
CCGTGGTGGGCGTTTCAAGCCCATCCCTGTTTGGGAGGCAAGGATCCCCAGTTCACTTCCGGCACTT
TCCAGATATATTTCAAGAGTGTGGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207235 protein sequence Red=Cloning site Green=Tags(s) MGPP LAPRPAHVPGEAGPRRTRESRPGAVSFADVAVYFSPEEWECLRPAQRALYRDMRETFGHLGALGF SVPKPAFISWVEGEVEAWSPEAQDPDGESSAAF SRGQGQEAGSRDGNEEKERLKKCPKQKEVAHEVAVKE WWPSVACPEFCNPRQSPMNPWLKDTLTRRLPHSCPDGGRNFSYPSLLASHQRVHSGERPFS CGQCQARFS QRRYLLQHQFIHTGEKPYPCDCGRRFRQGRSLAIHRRRAHTGEKPYACSDCKSRFTYPYLLAIHQKHTG EKPYS CPDCSLRFAYTSL LAIHRRRIHTGEKPYPCDCGRRFTYSLLLSHRRIHSDSRPFPCVECGKGFK RKTAL EAHRWIHRSCSERRAWQQA VVGRSEPIPV LGGKDP PVHFRHFPDIFQECG SGP TRTRPLEQKLI SEEDLA ANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152458
ORF Size:	1215 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_152458.4
RefSeq Size:	3265 bp
RefSeq ORF:	1218 bp
Locus ID:	146540
UniProt ID:	A8K8V0
Cytogenetics:	16p11.2
MW:	46.2 kDa