

Product datasheet for **RC207237L1V**

ZBPB2 (NM_198844) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZBPB2 (NM_198844) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZBPB2
Synonyms:	ZPBPL
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_198844
ORF Size:	948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207237).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198844.2
RefSeq Size:	2205 bp
RefSeq ORF:	951 bp
Locus ID:	124626
UniProt ID:	Q6X784
Cytogenetics:	17q21.1
Protein Families:	Secreted Protein
MW:	36.1 kDa



[View online »](#)

Gene Summary: Is implicated in sperm-oocyte interaction during fertilization.[UniProtKB/Swiss-Prot Function]