

Product datasheet for **RC207245L4V**

ILP2 (BIRC8) (NM_033341) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ILP2 (BIRC8) (NM_033341) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ILP2
Synonyms:	hILP2; ILP-2; ILP2; RNF136
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033341
ORF Size:	708 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207245).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033341.3 , NP_203127.3
RefSeq Size:	2047 bp
RefSeq ORF:	711 bp
Locus ID:	112401
UniProt ID:	Q96P09
Cytogenetics:	19q13.42
Protein Families:	Druggable Genome
MW:	26.9 kDa



[View online »](#)

Gene Summary: Protects against apoptosis mediated by BAX.[UniProtKB/Swiss-Prot Function]