

Product datasheet for **RC207247A1V**

Human ASB17 (NM_080868) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ASB17 (NM_080868) AAV Particle
Tag:	Myc-DDK
Symbol:	ASB17
Synonyms:	Asb-17
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207247 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTAAATCTACTAAATTATGTGGTAAGACTTCTTGTCGAAGAAGCAATATATTCTGCAATCTCCTTG
ACAAAATTGTTAAAAGACCCCTCCCTACAGTTTTTGGGTCAGTGGGGATATCACTGTTACGAACCAAGGAT
TTACAGATCACTGGCAAAAATCTGAGGTATGTGGACTTGGATGGTTTTGACGCACTACTCACAGATTAC
ATTGCATTTGTGAAAAATCAGGATACCGTTTTGAAGTAAGTTTAACTCGACTTCACTGAAATATGTG
TGAATACAATTCTGACTGGGTTTTGCCAGAAAAGGTAATCCTGACTTTGTGGAATTGCTTCTCAAGAA
GACAAAAGACTATGTTCAAGACAGAAGTTGTAACCTGGCACTGATATGGAGAAGTTTACACCAAGTATAC
TGTCGAAGCCCATTAAGTGCCATCACACCTCTCTTTATGTAGCTCAGACAAGACAGTCTAATATCTTCA
AAATACTACTGCAATATGGAATCTTAGAAAGAGAAAAAACCTATCAACATTGTCTTAACAATAGTACT
CTACCCTTCGAGAGTAAGAGTAATGGTTGATCGTGAATTGGCTGACATCCATGAAGATGCCAAAACATGT
TTGGTACTATGTTCCAGAGTGCTTTCTGTCATTTCAAGGAAATAAAGACACAGCTGAGTTTAGGAA
GACATCCAATTATTTCAAATTGGTTTGATTACATTCCTTCAACAAGATACAAAGATCCATGTGAAGTATT
ACATCTTTCAGACTAACCATCAGGAATCAACTATTAACCAACAATATGCTCCAGATGGAATATTTTCA
CTTCTAATTCCTGCTCGTCTACAAAACATCTGAATTTAGAAATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207247 protein sequence Red=Cloning site Green=Tags(s) MSKSTKLCGKTSCPRSNIFCNLDDKIVKRPSLQFLGQWGYHCYEPRIYRSLAKILRYVDLDGFDALLTDY IAFVEKSGYRFEVSFNLDFTICVNTILYWVFARKGNPDFVELLLKTKDYVQDRSCNLALIWRFTFPVY CPSPLSGITPLFYVAQTRQSNIFKILLQYGILEREKNPINIVLTIVLYPSRVVMVDRELADIHEDAKTC LVLCSSVLSVISVKEIKTQLSLGRHPPIISNWFYIPSTRYKDPCELLHLCRLTIRNQLLTNNMLPDGIFS LLIPARLQNYLNLEI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_080868
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080868.1</u>
RefSeq Size:	1121 bp
RefSeq ORF:	888 bp
Locus ID:	127247
UniProt ID:	<u>Q8WXJ9</u>
Cytogenetics:	1p31.1
MW:	34.3 kDa