

Product datasheet for **RC207247L4V**

ASB17 (NM_080868) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ASB17 (NM_080868) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ASB17
Synonyms:	Asb-17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080868
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080868.1
RefSeq Size:	1121 bp
RefSeq ORF:	888 bp
Locus ID:	127247
UniProt ID:	Q8WXJ9
Cytogenetics:	1p31.1
Protein Families:	Druggable Genome
MW:	34.3 kDa



[View online »](#)

Gene Summary:

May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]