

## Product datasheet for **RC207254L2V**

### HECTD2 (NM\_182765) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HECTD2 (NM_182765) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HECTD2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_182765
ORF Size:	2328 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207254).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182765.1</a>
RefSeq Size:	4963 bp
RefSeq ORF:	2331 bp
Locus ID:	143279
UniProt ID:	<a href="#">Q5U5R9</a>
Cytogenetics:	10q23.32
Protein Families:	Druggable Genome
MW:	88.1 kDa


[View online »](#)

**Gene Summary:**

Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.[UniProtKB/Swiss-Prot Function]