

## Product datasheet for **RC207257L4V**

### SPRN (NM\_001012508) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPRN (NM_001012508) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPRN
Synonyms:	bA108K14.1; SHADOO; SHO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001012508
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207257).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001012508.2</a> , <a href="#">NP_001012526.2</a>
RefSeq Size:	3188 bp
RefSeq ORF:	456 bp
Locus ID:	503542
UniProt ID:	<a href="#">Q5BIV9</a>
Cytogenetics:	10q26.3
MW:	14.5 kDa


[View online »](#)

**Gene Summary:**

Prion-like protein that has PrP(C)-like neuroprotective activity. May act as a modulator for the biological actions of normal and abnormal PrP (By similarity).[UniProtKB/Swiss-Prot Function]