

Product datasheet for **RC207258**

ZNF92 (NM_007139) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF92 (NM_007139) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF92
Synonyms:	HEL-203; HPF12; HTF12; TF12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC207258 ORF sequence, **codon optimized**.
 Due to the complexity of NM_007139, the ORF clone is codon optimized for mammalian Expression.
 The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**

ATGGTTGACAAGACCCAGTGATGTGCTCCCACTTCGCTCAGGACGTCTGGCCCGAACACAGCATTAAGGACTCTTTTCAGAAAGTCATACTCAGAACCTACGGGAAGTACGGGCATGAAAATCTTCAGCTTCGGAAAGACCAACAAGTCCGTGGATGCTTGAAGGTGTACAAGGCGGGTATAACGGGCTGAATCAGTGTTCGACTACA ACTGATCCAAGATTTTTCAGTGTGACAAATACGTTAAAGTGTTTCATAAATCCCAACGTTAATCGCAACAAGATCCGGCACACGGGGAAAAAACCTTCAAGTGTAAAGATCGCGGCAAAAGCTTTTGCATGTTGTCACAGCTCACTCAGCACAAAAAGATCCATACGCGAGAGTACAGTTACAAGTGTGAGGAATGTGGCAAGGCTTTAACTGGTCTAGCACCTTGACTAAACATAAGATTATACACACCGGGGAAAAACCTATAAGTGTGAGGAATGCGGAAAGGCTTTAATCGCTCTAGTAACCTCACAAAACACAAGATCATACACACCGGGGAGAAACCATACAAATGTGAAGAGTGTGGCAAGGCCTTTAACAGGTCTCCACTGACAAAGCATAAGCGGATTCACACTGAGGAGAAGCCTTACAAGTGCAGGAGTGCAGGCAAGGCTTTAACCAATTCTCTATCCTCAACAAACAAAAGGATCCATATGGAGGATAAACCTACAATGTGAAGAGTGTGGCAAAGCGTTTCGCGTCTTTCTATTCTGAAGAAACACAAAATCATCCACACCGGGGAAAAACCTATAAATGCGAAGAATGTGGCAAGGCA TTCAATCAGTTCAGTAACCTTGACCAAAACACAAGATCATTACACCGGGGAGAAAGCCATATAAGTGTGATG AATGCGGCAAGGCTTTCAACAGTCTAGCACGCTCACTAAGCATAAACGCATCCACACCGGCGAGAAGCCCTACAAATGCGAGGAGTGTGGCAAAGCTTTCAAGCAGTCCAGCACTCTGACCGAACATAAATCATCCACACCGGAGAAAAGCCTTATAAATGCGAAAAGTGCAGGCAAGGCGTTTTCTGGTCTTCAGCGTTTACTAAACATAAAAAGAAATCACATGGAAGACAAGCCATACAAGTGTGAAGAGTGCAGGAAAGCCTTTAGCGTCTTCTCTACCCTGACAAAACACAAGATCATCCACACGCGGAGAAGCCCTACAATGCGAGGAGTGCAGGCAAAAGCA TTTAATCAGTCTCTATCTTACAAAACACAAAATCATCCATACCGAAGGAAATCTTATAAATGTGAGAAGTGCAGCAATGCTTTCAACCAGAGCAGTAACCTGACTGCTAGGAAGATAATATACCGGAGAGAAACGTACAAGTACGAAGAGTGTGATAAAGCATTCAACAAGTTTTCAACGCTGATTACTACCAAAATCATTTACACCGGGGAAAAAGCCCTGTAACACGAATGCGGCGAGAGCATTCAATAAGAGCAGCAACTATACCAAGGAAA AACTGCAGACT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC207258 representing NM_007139
 Red=Cloning site Green=Tags(s)

MVDKTPVMCSHFAQDVWPEHSIKDSFQKVIILRTYGYGHENLQLRKDHKSVDACKVYKGGYNGLNQCLTTTDSKIFQCDKYVKVFHFKFPNVNRNKIRHTGKPKFKCKNRGKSFCMLSQLTQHKKIHTREYSYKCEECGKAFNWSSTLTKHKI IHTGEKPYKCEECGKAFNRSSNLTKHKI IHTGEKPYKCEECGKAFNRSSLTKHKRIHTEEKPYKCEECGKAFNQFSILNKHKRIHMEDKPYKCEECGKAFRVFSILKHKI IHTGEKPYKCEECGKAFNQFSNLTKHKI IHTGEKPYKCEECGKAFNQSSSTLTKHKRIHTGEKPYKCEECGKAFKQSSSTLTHEKIIHTGEKPYKCEKCGKAFSWSSAFTKHKRNHEDKPYKCEECGKAFSVFSTLTKHKI IHTREKPYKCEECGKAFNQSSIFTKHKI IHTEGSYKCEKCGNAFNQSSNLTARKIIYTGEKPYKYEEDKAFNKFSTLITHQIIYTGEKPKHECGRAFKNSSNYTKEKLQT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_007139

ORF Size: 1551 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007139.1](#), [NM_007139.2](#), [NM_007139.3](#), [NP_009070.2](#)

RefSeq Size: 3093 bp


RefSeq ORF: 1554 bp

Locus ID: 168374

UniProt ID: [Q03936](#)

Cytogenetics: 7q11.21
 Protein Families: Transcription Factors
 MW: 60.2 kDa
 Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC207258