

Product datasheet for **RC207259L3V**

RAB15 (NM_198686) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAB15 (NM_198686) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAB15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198686
ORF Size:	617 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207259).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198686.1
RefSeq Size:	3326 bp
RefSeq ORF:	627 bp
Locus ID:	376267
UniProt ID:	P59190
Cytogenetics:	14q23.3
MW:	23.3 kDa
Gene Summary:	May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.[UniProtKB/Swiss-Prot Function]



[View online »](#)