

Product datasheet for **RC207271**

SFRS5 (SRSF5) (NM_001039465) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SFRS5 (SRSF5) (NM_001039465) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SFRS5
Synonyms:	HRS; SFRS5; SRP40
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC207271 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTGGCTGTCGGGTATTCATCGGGAGACTAAATCCAGCGCCAGGGAGAAGGACGTGGAAAGATTCT
TCAAGGGATATGGACGGATAAGAGATATTGATCTGAAAAGAGGCTTTGGTTTTGTGGAATTTGAGGATCC
AAGGGATGCAGATGATGCTGTGTATGAGCTTGATGGAAAAGAAGCTCTGTAGTGAAAGGGTACTATTGAA
CATGCTAGGGCTCGTCCAGAGGTGGAAGAGGTAGAGGACGATACTCTGACCGTTTTAGTAGTCGCAGAC
CTCGAAATGATAGACGAAATGCTCCACCTGTAAGAACAGAAAATCGTCTTATAGTTGAGAATTTATCCTC
AAGAGTCAGCTGGCAGGATCTCAAAGATTTTCATGAGACAAGCTGGGGAAGTAACGTTTGGGATGCACAC
CGACCTAAATTAATGAAGGGTGGTTGAGTTGCCTCTTATGGTGACTTAAAGAATGCTATTGAAAAAC
TTTCTGAAAAGGAAATAAATGGGAGAAAAATAAATTAATTGAAGGCAGCAAAAGGCACAGTAGGTCAAG
AAGCAGGTCTCGATCCCGGACCAGAAGTTCCTTAGGTCTCGTAGCCGATCCCGTTCCCGTAGTCGCAAA
TCTTACAGCCGGTCAAGAAGCAGGAGCAGGAGCCGGAGCCGAGCAAGTCCCGTTCTGTTAGTAGGTCTC
CCGTGCCTGAGAAGGCCAGAAACGTGGTTCTTCAAGTAGATCTAAGTCTCCAGCATCTGTGGATCGCCA
GAGGTCCCGTCCCGATCAAGGTCCAGATCAGTTGACAGTGGCAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC207271 protein sequence
 Red=Cloning site Green=Tags(s)

MSGCRVFIGRLNPAAREKDVERFFKGYGRIRDIDLKRGFGVFEFEDPRDADDAVYELDGKELCSERVITIE
 HARARSRGGRGRGRYSDRFSSRRPRNDRRNAPPVRTENRLIVENLSSRVSWQDLKDFMRQAGEVTFADAH
 RPKLNEGVVAFASYGDLKNAIEKLSGKEINGRKIKLIEGSKRHSRSRSRSRTRSSRSRSRSRSRSRK
 SYSRSRSRSRSRSRSRSVSRSPVPEKSQKRGSSRSRSPASVDRQSRRSRSRSRSVDSGN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6377_b12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001039465

ORF Size: 816 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001039465.2
RefSeq Size:	1651 bp
RefSeq ORF:	819 bp
Locus ID:	6430
UniProt ID:	Q13243
Cytogenetics:	14q24.1
Protein Pathways:	Spliceosome
MW:	31.3 kDa
Gene Summary:	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]

Product images:



Circular map for RC207271



Western blot validation of overexpression lysate (Cat# [LY416317]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC218652] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).