

# Product datasheet for RC207276L3

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Monoamine Oxidase A (MAOA) (NM\_000240) Human Tagged Lenti ORF Clone

## **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** Monoamine Oxidase A (MAOA) (NM\_000240) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Monoamine Oxidase A

Synonyms: BRNRS; MAO-A

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-Myc-DDK-P2A-Puro (PS100092)

**E. coli Selection:** Chloramphenicol (34 ug/mL)

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC207276).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 






 $<sup>\</sup>ensuremath{^*}$  The last codon before the Stop codon of the ORF.





### Plasmid Map:



**ACCN:** NM\_000240 **ORF Size:** 1581 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 000240.2</u>

RefSeq Size: 5330 bp RefSeq ORF: 1584 bp



# Monoamine Oxidase A (MAOA) (NM\_000240) Human Tagged Lenti ORF Clone - RC207276L3

 Locus ID:
 4128

 UniProt ID:
 P21397

 Cytogenetics:
 Xp11.3

**Domains:** Amino\_oxidase

**Protein Families:** Druggable Genome

**Protein Pathways:** Arginine and proline metabolism, Drug metabolism - cytochrome P450, Glycine, serine and

threonine metabolism, Histidine metabolism, Metabolic pathways, Phenylalanine metabolism,

Tryptophan metabolism, Tyrosine metabolism

**MW:** 59.7 kDa

**Gene Summary:** This gene is one of two neighboring gene family members that encode mitochondrial

enzymes which catalyze the oxidative deamination of amines, such as dopamine,

norepinephrine, and serotonin. Mutation of this gene results in Brunner syndrome. This gene has also been associated with a variety of other psychiatric disorders, including antisocial behavior. Alternatively spliced transcript variants encoding multiple isoforms have been

observed. [provided by RefSeq, Jul 2012]