

## Product datasheet for **RC207281L3V**

### **PASD1 (NM\_173493) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PASD1 (NM_173493) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PASD1
Synonyms:	CT63; CT64; OXTES1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173493
ORF Size:	2319 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173493.1</a>
RefSeq Size:	3250 bp
RefSeq ORF:	2322 bp
Locus ID:	139135
UniProt ID:	<a href="#">Q8IV76</a>
Cytogenetics:	Xq28
MW:	87.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that is thought to function as a transcription factor. The protein is a cancer-associated antigen that can stimulate autologous T-cell responses, and it is therefore considered to be a potential immunotherapeutic target for the treatment of various hematopoietic malignancies, including diffuse large B-cell lymphoma. [provided by RefSeq, May 2010]