

Product datasheet for **RC207286L4V**

LDLRAD3 (NM_174902) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LDLRAD3 (NM_174902) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LDLRAD3
Synonyms:	LRAD3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_174902
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207286).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_174902.2
RefSeq Size:	2969 bp
RefSeq ORF:	1038 bp
Locus ID:	143458
UniProt ID:	Q86YD5
Cytogenetics:	11p13
Protein Families:	Druggable Genome, Transmembrane
MW:	37.4 kDa


[View online »](#)

Gene Summary:

May influence APP processing, resulting in a decrease in sAPP-alpha production and increased amyloidogenic P3 peptide production.[UniProtKB/Swiss-Prot Function]