

Product datasheet for RC207294L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MRCL3 (MYL12A) (NM_006471) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: MRCL3 (MYL12A) (NM_006471) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: MRCL3

Synonyms: HEL-S-24; MLC-2B; MLCB; MRCL3; MRLC3; MYL2B

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC207294).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_006471

ORF Size: 513 bp





MRCL3 (MYL12A) (NM_006471) Human Tagged Lenti ORF Clone - RC207294L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 006471.2</u>, <u>NP 006462.1</u>

 RefSeq Size:
 1243 bp

 RefSeq ORF:
 516 bp

 Locus ID:
 10627

 UniProt ID:
 P19105

 Cytogenetics:
 18p11.31

Domains: EFh

Protein Pathways: Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton, Tight

junction

MW: 19.8 kDa


Gene Summary: This gene encodes a nonsarcomeric myosin regulatory light chain. This protein is activated by

phosphorylation and regulates smooth muscle and non-muscle cell contraction. This protein may also be involved in DNA damage repair by sequestering the transcriptional regulator apoptosis-antagonizing transcription factor (AATF)/Che-1 which functions as a repressor of p53-driven apoptosis. Alternate splicing results in multiple transcript variants. A pseudogene


of this gene is found on chromosome 8.[provided by RefSeq, Dec 2014]



Product images:



Circular map for RC207294L1



Double digestion of RC207294L1 using Sgfl and Mlul