

Product datasheet for **RC207313L4V**

PLAC2 (TINCR) (NM_153375) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PLAC2 (TINCR) (NM_153375) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLAC2
Synonyms:	FLJ90734
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153375
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153375.1 , NP_699206.1
RefSeq Size:	2467 bp
RefSeq ORF:	455 bp
Locus ID:	257000
Cytogenetics:	19p13.3
MW:	16.5 kDa



[View online »](#)

Gene Summary:

This gene produces a spliced long non-coding RNA that binds RNAs. This transcript interacts with staufer-1 protein to regulate the stability of mRNAs for genes involved in the differentiation of epidermal tissue. Variation in this gene may be associated with cancer progression. [provided by RefSeq, Aug 2017]