

Product datasheet for **RC207316L1V**

MA2 (PNMA2) (NM_007257) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MA2 (PNMA2) (NM_007257) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MA2
Synonyms:	MA2; MM2; RGAG2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007257
ORF Size:	1092 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207316).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007257.4
RefSeq Size:	4788 bp
RefSeq ORF:	1095 bp
Locus ID:	10687
UniProt ID:	Q9UL42
Cytogenetics:	8p21.2
MW:	41.3 kDa


[View online »](#)