

Product datasheet for **RC207318L3V**

MMP15 (NM_002428) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MMP15 (NM_002428) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MMP15
Synonyms:	MMP-15; MT2-MMP; MT2MMP; MTMMP2; SMCP-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002428
ORF Size:	2007 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207318).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002428.2
RefSeq Size:	4456 bp
RefSeq ORF:	2010 bp
Locus ID:	4324
UniProt ID:	P51511
Cytogenetics:	16q21
Domains:	hemopexin, Peptidase_M10, ZnMc
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

MW: 75.8 kDa

Gene Summary: This gene encodes a member of the peptidase M10 family and membrane-type subfamily of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Members of this subfamily contain a transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. The encoded preproprotein is proteolytically processed to generate the mature protease. This protein may play a role in cancer progression. [provided by RefSeq, Jan 2016]