

Product datasheet for **RC207325A1V**

Human LPTM4B (NM_018407) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LPTM4B (NM_018407) AAV Particle
Tag:	Myc-DDK
Symbol:	LPTM4B
Synonyms:	LPTM4beta; LC27
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207325 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACGTACGGACTCGGGTCACATGGCCGAGTCCGCCCCGCCCTCCCGTCCCCGCCGCTGCAGCCG
TCGCCTTCGGAGCGAAGGATACCGACCCGGCAGAAGCTCGGAGCTCTCGGGTATCGAGGAGGCAGGCC
GCGGGCGCAGGGCGAGCGGGCCGGAGCCGGAGCGCGGAGGAGCCGGCAGCAGCGCGCGGGGCTC
CAGGCGAGGCGGTGACGCTCCTGAAAACCTGCGCGCGCGCTAGCGCCACTGCGCCGGAGCGATGAAGA
TGGTCGCGCCCTGGACGCGGTTCTACTCCAACAGCTGCTGCTTGTGCTGCCATGTCCGCACCGGCACCAT
CCTGCTCGGCGTCTGGTATCTGATCATCAATGCTGTGGTACTGTTGATTTTATTGAGTGCCTGGCTGAT
CCGGATCAGTATAACTTTTCAAGTTCTGAACTGGGAGGTGACTTTGAGTTCATGGATGATGCCAATATGT
GCATTGCCATTGCGATTTCTTCTCATGATCCTGATATGTGCTACGGCTACTTACGGAGCGTACAAGCA
ACGCGCAGCCTGGATCATCCCATTCTTCTGTTACCAGATCTTTGACTTTGCCCTGAACATGTTGTTGCA
ATCACTGTGCTTATTTATCCAACTCCATTGAGGAATACATACGGCAACTGCCTCCTAATTTCCCTACA
GAGATGATGTCATGTCAGTGAATCCTACCTGTTTGGTCCTTATTATTCTTCTGTTTATTAGCATTATCTT
GACTTTTAAGGGTTACTTGATTAGCTGTGTTTGAAGTCTACCGATACATCAATGGTAGGAACCTCTCT
GATGTCCTGGTTTATGTTACCAGCAATGACACTACGGTCTGCTACCCCGTATGATGATGCCACTGTGA
ATGGTCTGCCAAGGAGCCACCGCCACCTTACGTGTCTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207325 protein sequence Red=Cloning site Green=Tags(s) MTSRTRVTWSPPRPLVPAAAAVAFGAKGTDPAEARSSRGIEEAGPRAHGRAGREPERRRSRQQRGGL QARRSTLLKTCARASATAPGAMKMPWTRFYNSCCLCCHVRTGTILLGVWYLIINAVVLLILLSALAD PDQYNFSSSELGGDFEFMDANMCIAIAISLLMILICATATYGAYKQRAAWIIPFFCYQIFDFALNMLVA ITVLIYPNSIQEYIRQLPPNFPYRDDVMSVNPTCLVLIILLFISIILTFKGYLISCVWNCYRYINGRNSS DVLVYVTSNDTTVLLPPYDDATVNGAAKEPPPPYVSA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018407
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018407.4, NP_060877.3</u>
RefSeq Size:	2242 bp
RefSeq ORF:	681 bp
Locus ID:	55353
UniProt ID:	<u>Q86VI4</u>
Cytogenetics:	8q22.1
MW:	35 kDa