

## Product datasheet for **RC207341L2V**

### **C17ORF39 (GID4) (NM\_024052) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C17ORF39 (GID4) (NM_024052) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GID4
Synonyms:	C17orf39; VID2; VID24
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_024052
ORF Size:	900 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024052.4</a>
RefSeq Size:	4232 bp
RefSeq ORF:	903 bp
Locus ID:	79018
UniProt ID:	<a href="#">Q8IVV7</a>
Cytogenetics:	17p11.2
MW:	33.5 kDa



[View online »](#)

**Gene Summary:**

The multiprotein Mediator complex is a coactivator required for activation of RNA polymerase II transcription by DNA bound transcription factors. The protein encoded by this gene is thought to be a subunit of the Mediator complex. This gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]