

Product datasheet for **RC207345L3V**

P2RY5 (LPAR6) (NM_005767) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	P2RY5 (LPAR6) (NM_005767) Human Tagged ORF Clone Lentiviral Particle
Symbol:	P2RY5
Synonyms:	ARWH1; HYPT8; LAH3; LPA-6; P2RY5; P2Y5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005767
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207345).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005767.4
RefSeq Size:	3043 bp
RefSeq ORF:	1035 bp
Locus ID:	10161
UniProt ID:	P43657
Cytogenetics:	13q14.2
Domains:	7tm_1
Protein Families:	Druggable Genome, GPCR, Transmembrane


[View online »](#)

Protein Pathways: Neuroactive ligand-receptor interaction

MW: 39.4 kDa

Gene Summary: The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene aligns with an internal intron of the retinoblastoma susceptibility gene in the reverse orientation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009]