

## Product datasheet for **RC207368L2V**

### TMEM65 (NM\_194291) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM65 (NM_194291) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM65
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_194291
ORF Size:	720 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207368).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_194291.1</a> , <a href="#">NP_919267.1</a>
RefSeq Size:	3258 bp
RefSeq ORF:	723 bp
Locus ID:	157378
UniProt ID:	<a href="#">Q6PI78</a>
Cytogenetics:	8q24.13
Protein Families:	Transmembrane
MW:	25.3 kDa


[View online »](#)

**Gene Summary:**

May play an important role in cardiac development and function. May regulate cardiac conduction and the function of the gap junction protein GJA1. May contribute to the stability and proper localization of GJA1 to cardiac intercalated disk thereby regulating gap junction communication (By similarity). May also play a role in the regulation of mitochondrial respiration and mitochondrial DNA copy number maintenance (PubMed:28295037). [UniProtKB/Swiss-Prot Function]