

Product datasheet for **RC207401L4V**

ZNF398 (NM_020781) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF398 (NM_020781) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF398
Synonyms:	P51; P71; ZER6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020781
ORF Size:	1929 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207401).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020781.2
RefSeq Size:	5508 bp
RefSeq ORF:	1416 bp
Locus ID:	57541
UniProt ID:	Q8TD17
Cytogenetics:	7q36.1
Domains:	KRAB, zf-C2H2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 71.3 kDa

Gene Summary: This gene encodes a member of the Kruppel family of C2H2-type zinc-finger transcription factor proteins. The encoded protein acts as a transcriptional activator. Two transcript variants encoding distinct isoforms have been identified for this gene. Other transcript variants have been described, but their full length sequence has not been determined. [provided by RefSeq, Jul 2008]