

Product datasheet for **RC207441A1V**

Human ARHGEF39 (NM_032818) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ARHGEF39 (NM_032818) AAV Particle
Tag:	Myc-DDK
Symbol:	ARHGEF39
Synonyms:	C9orf100
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207441 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCTCTCCTGCCCCGTTTCGGGTGCCGGTGCAAGAGCAGCGTGCCCGCTGGGAGCGGAAACGCG
CCTGCACCGCCCGGGAGCTGCTAGAGACCGAGCGGCGCTACCAAGAACAGCTGGGGCTGGTGCCACGTA
CTTTTGGGGATCCTGAAAGCCAAGGGGACCCTGCGACCACCTGAGCGCCAGGCCCTGTTTGGCTCCTGG
GAGCTCATCTACGGCGCCAGCCAGGAGCTGCTCCCTACCTGGAAGGAGGATGCTGGGGCCAAGGGCTGG
AGGGCTTCTGCCGCCACTTGGAGCTCTATAACCAATTTGCTGCCAACTCAGAGAGGTCCCAGACCACCT
GCAGGAGCAGCTAAAGAAAAATAAAGGTTTCCGGAGGTTTGTACGGCTTCAGGAAGGCCGCCCTGAGTTT
GGGGGCCCTTCAGCTCCAGGACCTGCTCCCTCTGCCTCTGCAACGGCTCCAGCAGTATGAGAATCTCGTCG
TAGCTTTGGCTGAAAACACAGGTCCCAACAGCCCTGACCATCAACAGCTCACACGGGCTGCCGACTGAT
AAGTGAGACTGCCAGAGAGTCCATACTATTGGTCAGAAACAGAAGAATGACCAGCACCTTCGGCGTGTC
CAGGCTCTGCTCAGTGGACGCCAGGCAAGGGGCTGACCTCAGGGCGCTGGTTCCTACGCCAGGGCTGGC
TGTTAGTGGTGCCTCCCATGGGAGCCTCGGCCCGCATGTTCTTCTTCACTGATGTGCTCCTCAT
GGCCAAGCCTCGGCCTCACTGCACCTGCTGCGGAGTGGCACCTTTGCCTGCAAGGCCCTTACCCCATG
GCCAGTGTCTCAGCAGGGTCTTTGGCCACTCAGGAGGCCCTTGTTGGTGGTTGCTCAGTCTGTCT
TCCCTCATGAGAAGCTACTGCTTATGTCCACAGACCAGGAGGAGCTGTACGCTGGTACCACAGTCTGAC
TTGGGCTATCAGCAGCCAGAAAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207441 protein sequence Red=Cloning site Green=Tags(s) MELSCPGSRCPVQEQRARWERKRACTARELLETERRYQEQLGLVATYFLGILKAKGTLRPPPERQALFGSW ELIYGASQELLPYLEGGCWGQGLEGFGRHLELYNQFAANSERSQTTLQEQLKKNKGFRFVRLQEGRPEF GGLQLQDLLPLPLQRLQQYENLVVALAENTGPNSPDHQQLTAAARLISETAQRVHTIGQKQKNDQHLRRV QALLSGRQAKGLTSGRWFLRQGWLLVPPPHGEPRPMFFLFTDVLLMAKPRPPLHLLRSGTFACKALYPM AQCHLSRVFGHSGGPCGGLLSLSPHEKLLLMSTDQEELSRWYHSLTWAISSQKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032818
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032818.2</u>
RefSeq Size:	3763 bp
RefSeq ORF:	1008 bp
Locus ID:	84904
UniProt ID:	<u>Q8N4T4</u>
Cytogenetics:	9p13.3
MW:	38.3 kDa