

Product datasheet for **RC207441L3V**

ARHGEF39 (NM_032818) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARHGEF39 (NM_032818) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGEF39
Synonyms:	C9orf100
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032818
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207441).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032818.2
RefSeq Size:	3763 bp
RefSeq ORF:	1008 bp
Locus ID:	84904
UniProt ID:	Q8N4T4
Cytogenetics:	9p13.3
MW:	38.3 kDa
Gene Summary:	Promotes cell proliferation.[UniProtKB/Swiss-Prot Function]


[View online »](#)