

Product datasheet for **RC207445**

RPS15A (NM_001030009) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPS15A (NM_001030009) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RPS15A
Synonyms: DBA20; S15a
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC207445 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGCGCATGAATGTCCTGGCAGATGCTCTCAAGAGTATCAACAATGCCGAAAAGAGAGGCAAACGCC
AGGTGCTTATTAGGCCGTGCTCCAAAGTCATCGTCCGTTTCTCACTGTGATGATGAAGCATGGTTACAT
TGGCGAATTTGAAATCATTGATGACCACAGAGCTGGGAAAATTGTTGTGAACCTCACAGGCAGGCTAAAC
AAGTGTGGGTGATCAGCCCAGATTTGACGTGCAACTCAAAGACCTGGAAAAATGGCAGAATAATCTGC
TTCCATCCGCCAGTTTGGTTTCATTGTACTGACAACCTCAGCTGGCATCATGGACCATGAAGAAGCAAG
ACGAAAACACACAGGAGGAAAATCCTGGGATTCTTTTTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207445 protein sequence
Red=Cloning site Green=Tags(s)

MVRMNVLADALKSINNAEKRGKRQVLIRPCSKVIVRFLTYMMKHGYIGEFEIIDDHRAGKIVVNLTGRLN
KCGVISPRFDVQLKDLEKWNLLPSRQFGFIVLTTTSAGIMDHEEARRKHTGGKILGFFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6021_g03.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001030009

ORF Size: 390 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001030009.2](#)

RefSeq Size: 488 bp

RefSeq ORF: 393 bp


Locus ID: 6210

UniProt ID: [P62244](#)


Cytogenetics: 16p12.3
Protein Pathways: Ribosome
MW: 14.8 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S8P family of ribosomal proteins. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC207445



Western blot validation of overexpression lysate (Cat# [LY400408]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC207445 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).