

Product datasheet for **RC207446**

LSMD1 (NAA38) (NM_032356) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: LSMD1 (NAA38) (NM_032356) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: LSMD1
Synonyms: LSMD1; MAK31; PFAAP2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC207446 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGTTGCAGTCGGCGTCAGAGCAGCTCCAGTGCTGGGGTTAGCACGGGCTCTGGTTCTGGGACTGA
GGGCGAGTCAAGCAGCAAGATGGAGGGGTGGGAACTGCAGCTCCCGGCAGCCTCTGGGCTTTGTGTGA
GTGCCCGCGGGTCCCGGGAGCTGTGGTTCGCGGGCGGGCAGCCTGGGCTGCGCGGTCTGAGCGCTTG
TTACACCCCGCGCAGGATTCGGACGGAGAGCGCAGGACTCGGCGGTGAGCGCGCCGACAGCAGCTAG
AGGCGCTGCTCAACAAGACTATGCGCATTGCGATGACAGATGGACGGACTGGTGGCTGCTTCTCTG
CACTGACCGTGACTGCAATGTATCCTGGGCTCGGCGCAGGAGTTCCTCAAGCCGTGGATTCTTCTCT
GCCGGGAGCCCCGTGTGCTGGGCTGGCCATGGTACCCGGACACCACATCGTTTCCATTGAGGTGCAGA
GGGAGAGTCTGACCGGGCCTCCGTATCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207446 protein sequence
Red=Cloning site Green=Tags(s)

MAVAVGVRAAPVLGLARALVLGLRGSQAARWRWGTAAPGSLWALCECPAGCRELWFRGRAAWAARSERL
LHPAQSDGEREDSAAERARQQLLEALLNKTMIRMTDGRTL VGCF LCTDRDCNVILGSAQEF LKPSDFS
AGEPRVLGLAMVPGHHIVSIEVQRESLTGPPYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6335_g09.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_032356

ORF Size: 519 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032356.5](#)

RefSeq Size: 1023 bp

RefSeq ORF: 522 bp

Locus ID: 84316


UniProt ID: [Q9BRA0](#)

Cytogenetics: 17p13.1

MW: 18.7 kDa

Gene Summary: Auxillary component of the N-terminal acetyltransferase C (NatC) complex which catalyzes acetylation of N-terminal methionine residues.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC207446