

Product datasheet for **RC207456L2V**

RASSF4 (NM_032023) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RASSF4 (NM_032023) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RASSF4
Synonyms:	AD037
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_032023
ORF Size:	963 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207456).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032023.3
RefSeq Size:	2509 bp
RefSeq ORF:	966 bp
Locus ID:	83937
UniProt ID:	Q9H2L5
Cytogenetics:	10q11.21
Domains:	RA
Protein Families:	Druggable Genome



[View online »](#)

MW: 36.7 kDa

Gene Summary: The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. [provided by RefSeq, Jul 2008]