

Product datasheet for **RC207479L4V**

OTX2 (NM_021728) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	OTX2 (NM_021728) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OTX2
Synonyms:	CPHD6; MCOP55
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021728
ORF Size:	891 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207479).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021728.2
RefSeq Size:	2219 bp
RefSeq ORF:	894 bp
Locus ID:	5015
UniProt ID:	P32243
Cytogenetics:	14q22.3
Domains:	homeobox, TF_Otx



[View online »](#)

Protein Families:	Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency, Transcription Factors
MW:	32.4 kDa
Gene Summary:	This gene encodes a member of the bicoid subfamily of homeodomain-containing transcription factors. The encoded protein acts as a transcription factor and plays a role in brain, craniofacial, and sensory organ development. The encoded protein also influences the proliferation and differentiation of dopaminergic neuronal progenitor cells during mitosis. Mutations in this gene cause syndromic microphthalmia 5 (MCOPS5) and combined pituitary hormone deficiency 6 (CPHD6). This gene is also suspected of having an oncogenic role in medulloblastoma. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Pseudogenes of this gene are known to exist on chromosomes two and nine. [provided by RefSeq, Jul 2012]

Product images:



[RC207479L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC207479L4V particle to overexpress human OTX2-mGFP fusion protein.