

Product datasheet for **RC207480L4V**

TXNDC13 (TMX4) (NM_021156) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TXNDC13 (TMX4) (NM_021156) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TXNDC13
Synonyms:	DJ971N18.2; PDIA14; TXNDC13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021156
ORF Size:	1047 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207480).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021156.2 , NP_066979.2
RefSeq Size:	2368 bp
RefSeq ORF:	1050 bp
Locus ID:	56255
UniProt ID:	Q9H1E5
Cytogenetics:	20p12.3
Protein Families:	Druggable Genome, Transmembrane
MW:	38.8 kDa


[View online »](#)

Gene Summary:

This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, a catalytically active thioredoxin domain, one transmembrane domain and C-terminal ASP/GLU-rich calcium binding domain. Unlike most members of this gene family, it lacks a C-terminal ER-retention sequence. The encoded protein has been shown to have reductase activity in vitro. [provided by RefSeq, Jan 2017]