

## Product datasheet for **RC207485L2V**

### Scramblase 1 (PLSCR1) (NM\_021105) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Scramblase 1 (PLSCR1) (NM_021105) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Scramblase 1
Synonyms:	MMTRA1B
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_021105
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021105.1</a>
RefSeq Size:	2228 bp
RefSeq ORF:	957 bp
Locus ID:	5359
UniProt ID:	<a href="#">O15162</a>
Cytogenetics:	3q24
Domains:	Scramblase
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)

**MW:** 35 kDa

**Gene Summary:** May mediate accelerated ATP-independent bidirectional transbilayer migration of phospholipids upon binding calcium ions that results in a loss of phospholipid asymmetry in the plasma membrane. May play a central role in the initiation of fibrin clot formation, in the activation of mast cells and in the recognition of apoptotic and injured cells by the reticuloendothelial system.[UniProtKB/Swiss-Prot Function]