

Product datasheet for **RC207500L3V**

TMEM126B (NM_018480) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM126B (NM_018480) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM126B
Synonyms:	HT007; MC1DN29
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018480
ORF Size:	600 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207500).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018480.2 , NP_060950.2
RefSeq Size:	1038 bp
RefSeq ORF:	693 bp
Locus ID:	55863
UniProt ID:	Q8IUX1
Cytogenetics:	11q14.1
Protein Families:	Transmembrane
MW:	22.8 kDa



[View online »](#)

Gene Summary:

This gene encodes a mitochondrial transmembrane protein which is a component of the mitochondrial complex I assembly complex. The encoded protein serves as an assembly factor that is required for formation of the membrane arm of the complex. It interacts with NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 13. Naturally occurring mutations in this gene are associated with isolated complex I deficiency. A pseudogene of this gene has been defined on chromosome 9. [provided by RefSeq, Apr 2017]