

Product datasheet for **RC207509A1V**

Human Keratocan (KERA) (NM_007035) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Keratocan (KERA) (NM_007035) AAV Particle
Tag:	Myc-DDK
Symbol:	Keratocan
Synonyms:	CNA2; KTN; SLRR2B
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC207509 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGGCACAATCTGTTTCATCATGTGGGTGTTATTCATAACAGACACTGTGTGGTCTAGAAGTGTAAGGCAGGTCTATGAAGTACATGATTGAGTGGACTATTCATGACTTCGAGTGTCCCATGGAATGTTCTGCCCCACCCAGTTTCTACTGCTTTATTTGTGAAATAGAGGTCTCAAAGAAATTCCTGCTATTCTCTCAAGAATTTGGTATCTTTATCTTCAAAACAACCTGATAGAAACCATTCTGAAAAGCCATTTGAGAATGCCACCCAGCTAAGATGGATAAATCTAAACAAGAACAAATAACCAACTACGGAATTGAAAAGGAGCCCTAAGCCAGCTGAAGAAGTTGCTCTTCTTATTTCTGGAAGATAATGAGCTAGAGGAGGTACCTTCTCCATTGCCAAGAAGTTTGAACAATTACAATTAGCTAGAAATAAGGTGTCCAGAATTCCTCAAGGGACCTTTAGCAATCTGGAGAACCTGACCCTTCTTGACCTACAGAACAAATAGTGGACAATGCCTTTCAAAGAGACACTTTAAAGGACTCAAGAATCTCATGCAGCTAACATGGCCAAGAATGCCCTGAGGAATATGCCTCCAAGATTACCAGCCAATACAATGCAGTTGTTTTAGACAACAATTCATTGAAGGAATACCAGAAAATTATTTAATGTGATTCTAAAGTGGCCTTTTGGAGCTAAATCACAACAACTGTCAGATGAGGGTCTCCCATCAAGAGGATTTGATGTATCATCAATTCTAGATCTTCAACTGTCGCACAATCACTCACAAGGTTCCCGAATCAGTGCTCATCTGCAGCACCTTACCTTGATCATAACAAATTAAGTGTGAATGTCTCTGTAATATGTCACAGCCATCCATGCTGCCTGCAGAACGAGATTCTTCAGTTATGGACCTCATCTCGCTACCTCCGTCTGGATGGAATGAAATCAAACCACCAATTCGAATGGCTTTAATGACCTGCTTCAGACTTCTGCAGGCTGTCATTATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC207509 protein sequence Red=Cloning site Green=Tags(s) MAGTICFIMWVLFITDTVWSRSVRQVVEVHDSDDWTIHDFECPMECFPPSFPTALYCENRGLKEIPAIP SRIWYLYLQNNLIETIPEKPFENATQLRWINLNKNKITNYGIEKGALSQLKLLFLFLEDNEEEVPSPL PRSLEQLQLARNKVSRIPOGTFSNLENLTLLDLQNNKLVDNAFQRDTFKGLKNLMQLNMAKNALRNMPPR LPANTMQLFLDNNNSIEGIPENYFNVIPKVAFLRLNHNKLSDEGLPSRGFDVSSILDQLSHNQLTKVPRI SAHLQHLHLDHNKIKSVNVSVICSPSMLPAERDSFSYGPHRLRYLRDGNIEKPPIPMALMTCFRLLQAV II TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_007035
ORF Size:	1056 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_007035.3
RefSeq Size:	2543 bp
RefSeq ORF:	1059 bp
Locus ID:	11081
UniProt ID:	O60938
Cytogenetics:	12q21.33
MW:	40.5 kDa