

Product datasheet for **RC207513L4V**

RBMX2 (NM_016024) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RBMX2 (NM_016024) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBMX2
Synonyms:	CGI-79; Snu17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016024
ORF Size:	966 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207513).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016024.2
RefSeq Size:	1528 bp
RefSeq ORF:	969 bp
Locus ID:	51634
UniProt ID:	Q9Y388
Cytogenetics:	Xq26.1
MW:	37.3 kDa



[View online »](#)

Gene Summary: Involved in pre-mRNA splicing as component of the activated spliceosome.[UniProtKB/Swiss-Prot Function]