

Product datasheet for **RC207516L4V**

BET1L (NM_016526) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BET1L (NM_016526) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BET1L
Synonyms:	BET1L1; GOLIM3; GS15; HSPC197
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016526
ORF Size:	306 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207516).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016526.3
RefSeq Size:	2971 bp
RefSeq ORF:	309 bp
Locus ID:	51272
Cytogenetics:	11p15.5
Protein Families:	Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport
MW:	11.1 kDa



[View online »](#)

Gene Summary:

Vesicle SNARE required for targeting and fusion of retrograde transport vesicles with the Golgi complex. Required for the integrity of the Golgi complex (By similarity).
[UniProtKB/Swiss-Prot Function]