

## Product datasheet for **RC207524L2V**

### Neuroigin 1 (NLGN1) (NM\_014932) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Neuroigin 1 (NLGN1) (NM_014932) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Neuroigin 1
Synonyms:	NL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_014932
ORF Size:	2469 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207524).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014932.2</a>
RefSeq Size:	4947 bp
RefSeq ORF:	2472 bp
Locus ID:	22871
UniProt ID:	<a href="#">Q8N2Q7</a>
Cytogenetics:	3q26.31
Domains:	COesterase
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Protein Pathways:** Cell adhesion molecules (CAMs)

**MW:** 92 kDa

**Gene Summary:** This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. [provided by RefSeq, Jul 2008]