

## Product datasheet for **RC207526L4V**

### SSX2IP (NM\_014021) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SSX2IP (NM_014021) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SSX2IP
Synonyms:	ADIP; hMsd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014021
ORF Size:	1842 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207526).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014021.2</a>
RefSeq Size:	5658 bp
RefSeq ORF:	1845 bp
Locus ID:	117178
UniProt ID:	<a href="#">Q9Y2D8</a>
Cytogenetics:	1p22.3
Protein Pathways:	Adherens junction
MW:	71.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that binds the cancer-testis antigen Synovial Sarcoma X breakpoint 2 protein. The encoded protein may regulate the activity of Synovial Sarcoma X breakpoint 2 protein in malignant cells. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 3. [provided by RefSeq, Oct 2009]