

## Product datasheet for **RC207582L4V**

### RGL2 (NM\_004761) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	RGL2 (NM_004761) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RGL2
Synonyms:	HKE1.5; KE1.5; RAB2L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004761
ORF Size:	2331 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207582).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004761.2</a>
RefSeq Size:	3362 bp
RefSeq ORF:	2334 bp
Locus ID:	5863
UniProt ID:	<a href="#">O15211</a>
Cytogenetics:	6p21.32
Domains:	RA, RasGEFN, RasGEF
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 83.5 kDa

**Gene Summary:** Probable guanine nucleotide exchange factor. Putative effector of Ras and/or Rap. Associates with the GTP-bound form of Rap 1A and H-Ras in vitro (By similarity).[UniProtKB/Swiss-Prot Function]