

Product datasheet for **RC207583L3V**

SPARCL1 (NM_004684) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SPARCL1 (NM_004684) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPARCL1
Synonyms:	MAST 9; MAST9; PIG33; SC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004684
ORF Size:	1992 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207583).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004684.2
RefSeq Size:	2910 bp
RefSeq ORF:	1995 bp
Locus ID:	8404
UniProt ID:	Q14515
Cytogenetics:	4q22.1
Domains:	kazal, FOLN
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

MW: 75.2 kDa