

Product datasheet for **RC207586L4V**

CIAO1 (NM_004804) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CIAO1 (NM_004804) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CIAO1
Synonyms:	CIA1; WDR39
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004804
ORF Size:	1017 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207586).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004804.2
RefSeq Size:	4053 bp
RefSeq ORF:	1020 bp
Locus ID:	9391
UniProt ID:	O76071
Cytogenetics:	2q11.2
Domains:	WD40
Protein Families:	Transcription Factors



[View online »](#)

MW: 37.8 kDa

Gene Summary: Key component of the cytosolic iron-sulfur protein assembly (CIA) complex, a multiprotein complex that mediates the incorporation of iron-sulfur cluster into extramitochondrial Fe/S proteins (PubMed:17937914, PubMed:23891004). As a CIA complex component, interacts specifically with CIAO2A or CIAO2B and MMS19 to assist different branches of iron-sulfur protein assembly, depending of its interactors. The complex CIAO1:CIAO2B:MMS19 binds to and facilitates the assembly of most cytosolic-nuclear Fe/S proteins. CIAO1:CIAO2A specifically matures ACO1 and stabilizes IREB2 (PubMed:23891004). Seems to specifically modulate the transactivation activity of WT1 (PubMed:9556563). As part of the mitotic spindle-associated MMXD complex it may play a role in chromosome segregation (PubMed:20797633). [UniProtKB/Swiss-Prot Function]