

Product datasheet for **RC207599A1V**

Human TFF1 (NM_003225) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TFF1 (NM_003225) AAV Particle
Tag: Myc-DDK
Symbol: TFF1
Synonyms: BCEI; D21S21; HP1.A; HPS2; pNR-2; pS2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC207599 representing NM_003225
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCACCATGGAGAACAAGGTGATCTGCGCCCTGGTCTGGTGTCCATGCTGGCCCTCGGCACCCTGG
CCGAGGCCAGACAGAGACGTGTACAGTGGCCCCCGTGAAAGACAGAATTGTGGTTTCCTGGTGTAC
GCCCTCCAGTGTGCAAATAAGGGCTGCTGTTTCGACGACACCGTTCTGTGGGTCCCCTGGTCTCTAT
CCTAATACCATCGACGTCCCTCCAGAAGAGGAGTGTGAATTT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC207599 representing NM_003225
Red=Cloning site Green=Tags(s)

MATMENKVICALVLVSMALGLTAEQTTCTVAPRRQNCGFPGVTPSQCANKGCCFDDTVRGVPWCFY
PNTIDVPPEEECEF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_003225
ORF Size: 252 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003225.2</u>
RefSeq Size:	508 bp
RefSeq ORF:	255 bp
Locus ID:	7031
UniProt ID:	<u>P04155</u>
Cytogenetics:	21q22.3
MW:	9.15 kDa