

Product datasheet for **RC207622L4V**

D aspartate oxidase (DDO) (NM_003649) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	D aspartate oxidase (DDO) (NM_003649) Human Tagged ORF Clone Lentiviral Particle
Symbol:	D aspartate oxidase
Synonyms:	DASOX; DDO-1; DDO-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003649
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207622).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003649.2
RefSeq Size:	1733 bp
RefSeq ORF:	1110 bp
Locus ID:	8528
UniProt ID:	Q99489
Cytogenetics:	6q21
Domains:	DAO
Protein Pathways:	Alanine, aspartate and glutamate metabolism



[View online »](#)

MW: 41 kDa

Gene Summary: The protein encoded by this gene is a peroxisomal flavoprotein that catalyzes the oxidative deamination of D-aspartate and N-methyl D-aspartate. Flavin adenine dinucleotide or 6-hydroxyflavin adenine dinucleotide can serve as the cofactor in this reaction. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2019]