

## Product datasheet for **RC207641L4V**

### TCP10L (NM\_144659) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TCP10L (NM_144659) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TCP10L
Synonyms:	C21orf77; PRED77; TCP10A-1; TCP10A-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144659
ORF Size:	645 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207641).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_144659.3</a> , <a href="#">NP_653260.1</a>
RefSeq Size:	2730 bp
RefSeq ORF:	648 bp
Locus ID:	140290
UniProt ID:	<a href="#">Q8TDR4</a>
Cytogenetics:	21q22.11
MW:	23.9 kDa



[View online »](#)

**Gene Summary:**

May be involved in transcriptional regulation. Has in vitro transcription inhibition activity. Acts as a tumor suppressor in hepatocellular carcinoma (HCC) cells.[UniProtKB/Swiss-Prot Function]