

Product datasheet for **RC207668L2V**

FUK (FCSK) (NM_145059) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FUK (FCSK) (NM_145059) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FCSK
Synonyms:	1110046B12Rik; CDGF2; FUK
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_145059
ORF Size:	3252 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207668).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145059.2
RefSeq Size:	3923 bp
RefSeq ORF:	3255 bp
Locus ID:	197258
UniProt ID:	Q8N0W3
Cytogenetics:	16q22.1
Domains:	GHMP_kinases
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways

MW: 117.6 kDa

Gene Summary: The protein encoded by this gene belongs to the GHMP (galacto-, homoserine, mevalonate and phosphomevalonate) kinase family and catalyzes the phosphorylation of L-fucose to form beta-L-fucose 1-phosphate. This enzyme catalyzes the first step in the utilization of free L-fucose in glycoprotein and glycolipid synthesis. L-fucose may be important in mediating a number of cell-cell interactions such as blood group antigen recognition, inflammation, and metastasis. While several transcript variants may exist for this gene, the full-length nature of only one has been described to date. [provided by RefSeq, Jul 2008]