

## Product datasheet for **RC207672L2V**

### IL13 receptor alpha 2 (IL13RA2) (NM\_000640) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL13 receptor alpha 2 (IL13RA2) (NM_000640) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IL13 receptor alpha 2
Synonyms:	CD213A2; CT19; IL-13R; IL13BP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000640
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC207672).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000640.2</a>
RefSeq Size:	1376 bp
RefSeq ORF:	1143 bp
Locus ID:	3598
UniProt ID:	<a href="#">Q14627</a>
Cytogenetics:	Xq23
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Jak-STAT signaling pathway



[View online »](#)

**MW:** 44.2 kDa

**Gene Summary:** The protein encoded by this gene is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13. [provided by RefSeq, Jul 2008]